



Medela Swing Single and Medela Swing Maxi Diagram



This diagram shows what is in each bag and images of parts for visual reference. For more information refer to the Medela Pump instruction booklets on Medela.com

Medal Swing Single	Medela Swing Double
<p>1 x PersonalFit Breast Shield size M (24 mm) 1 x 2-component connector (1 x valve head and 1 x valve membrane) 1 x bottle 1 x bottle stand 1 x Calma Solitaire 1 x Calma Multi lid 1 x Swing tubing 1 x motor 1 x Main adaptor 1 x Swing storage bag 1 x sling (to hang pump from neck)</p>  <p>1x PersonalFit breast shield M (24mm) (for other breast shield sizes see Chapter 12, page 28)</p> <p>1x 2-component connector</p> <p>1x valve head</p> <p>2x valve membrane</p> <p>1x bottle</p> <p>1x motor unit</p> <p>On/Off button</p> <p>LED</p> <p>Reduce vacuum</p> <p>Increase vacuum</p> <p>Let-down button</p> <p>1x Swing tubing</p> <p>1 x bottle stand</p> <p>1x Calma Solitaire</p> <p>1x Calma multi lid</p> <p>1x mains adapter</p> <p>1x sling</p>	<p>2 x PersonalFit Breast shields size M(24 mm) 2 x Freestyle connectors (2 x breast shield bases, 2 x membranes, 2x transparent caps) 2 x bottles 2 x bottle stands 2 x Calma multi lid 1 x Calma Solitaire 1 x Motor unit 1 x Swing Maxi Tubbing 1 x Main Adaptor</p>  <p>2x PersonalFit™ Breastshield M (24mm) (for other breast shield sizes see Chapter 12, page 30)</p> <p>2x Freestyle connectors (breast shield base, membrane, transparent cap)</p> <p>2x bottle</p> <p>2x bottle stand</p> <p>1x Calma Solitaire</p> <p>2 x Calma Multi lid</p> <p>1x mains adapter</p> <p>1x motor unit</p> <p>On/Off button</p> <p>LED</p> <p>Reduce vacuum</p> <p>Increase vacuum</p> <p>Let-down button</p> <p>Short tubing</p> <p>Tubing support</p> <p>Long tubing</p> <p>1x Swing Maxi tubing</p>